

Access statement: equipment and measurements

A single-storey adapted bungalow in Whitstable, Kent, for disabled adults, families and carers. This is the spec sheet for occupational therapists, case managers and funding panels: makes, models, safe working loads and clearances, taken from the manufacturers' own spec sheets and the door measurements we have published.

965 mm

Front door, clear opening

926 mm

Internal doors, clear width

180 kg

Ceiling track hoist, safe working load

IN THIS STATEMENT [Getting in](#) [Beds and hoists](#) [Wet room](#) [Kitchen](#) [What we can't promise](#)

Getting in

The driveway is resin-bound, level, and takes two cars including an adapted vehicle. There is no step up into the house from where you park.

Parking and getting around

Single storey, no internal steps

DRIVEWAY

Private, resin-bound: two off-road spaces, adapted vehicles welcome

STREET PARKING

Outside the house if you need extra room. No residents permit on this road; check signs on arrival in case street rules change

CAR TO FRONT DOOR

Step-free path

FLOORS

Level flooring throughout, no carpet lips

PATIO

Level hard-standing suitable for wheelchairs, with the enclosed garden and BBQ area beyond

Doorways and thresholds		965 mm
Clear openings, as published on our doorway graphic		Front door clear opening
FRONT DOOR	965 mm clear, level threshold, no step	
PORCH OUTER DOORS	1720 mm opening, with a full-width ribbed entrance mat	
ALL INTERNAL DOORS	926 mm clear	
CONSERVATORY FRENCH DOORS	1720 mm opening	
REAR FRENCH DOORS	1720 mm opening, onto the patio	
PATIO THRESHOLD RAMP	20 mm high × 80 mm wide, non-slip	
DOOR FURNITURE	White doors with black lever handles, for contrast	
PORTABLE RAMPS	Fold-up ramps kept for the front door, yours for the stay and for days out	

Beds and hoists

The accessible bedroom can be set up with one or two profiling beds, depending on what you need.

Please bring your own slings

Slings are prescribed items of personal care. They are fitted to the individual and to the hoist mechanism, so we cannot safely supply them. Both hoists here take loop-style slings on a two-point spreader bar. Bring the sling the person already uses.

Amico GoLift 400: ceiling track hoist		180 kg
Main bedroom, over the profiling bed		Safe working load
SAFE WORKING LOAD	180 kg / 28 st	
LIFT AND LOWER	Powered	
TRAVERSE	Manual	
TROLLEY	Fixed	
ATTACHMENT	Lifting strap with loop	

Amico GoLift 400: ceiling track hoist Main bedroom, over the profiling bed		180 kg Safe working load
COVERAGE	Full-room	
THOROUGH EXAMINATION	LOLER, every six months	

Accora CommunityBed Accessible bedroom, under the ceiling track		180 kg Maximum user weight
SAFE WORKING LOAD	215 kg / 34 st	
MAXIMUM USER WEIGHT	180 kg / 28 st	
MATTRESS PLATFORM	900 × 2000 mm	
OVERALL FRAME	930 × 2250 mm	
PLATFORM HEIGHT	220–800 mm	
PROFILING	Four-section: backrest with auto-regression, kneebreak, Trendelenburg and anti-Trendelenburg	
MATTRESS	Accora Allevia Comfort FirmEdge 90 cm, or Roma Medical MATT 1, both pressure-relieving	
SIDE RAILS	Fabric siderail set, 200 cm	
MINIMUM USER HEIGHT	1460 mm	
CARER ACCESS	Space for a carer on both sides of the bed	

Joerns Oxford Midi 180: mobile hoist Kept in the bungalow, usable in any room		180 kg Safe working load
SAFE WORKING LOAD	180 kg / 28 st	
TURNING RADIUS	1235 mm	
LEGS OPEN, EXTERNAL	1170 mm	
LEGS CLOSED, EXTERNAL	600 mm	
LEG HEIGHT	100 mm	
GROUND CLEARANCE	25 mm	
SPREADER BAR HEIGHT	525–1660 mm	

Joerns Oxford Midi 180: mobile hoist

Kept in the bungalow, usable in any room

180 kg

Safe working load

EMERGENCY DESCENT

Manual

THOROUGH EXAMINATION

LOLER, every six months. Last examined 14 May 2026

Bringing your own mobile hoist

The bed platform lowers to 220 mm. The Oxford Midi's legs stand 100 mm high with 25 mm of ground clearance, and its turning radius is 1235 mm. If you are bringing your own hoist, measure its leg height and turning circle against those figures, and tell us what you are bringing so we can set the room up for it.

AAL RS4: transfer and standing aid

185 kg

Safe working load

The RS4 is an AAL stand assist. It is not an Arjo Sara Stedy, and the two are not interchangeable in a handling plan. Please check your care plan names the right one.

Wet room

Level-access throughout, with no lip at the entrance. Specified by Care Spaces.

Level-access wet room

Specified by Care Spaces

FLOOR

Level access with no lip, floor-level drain

TOILET SUPPORT RAIL

NYMAS hinged lift-and-lock, with drop-down leg

ASSISTANCE ALARM

Pull-cord with reset plate and over-door light

EXTRACTOR FAN

Fitted

RAZ-AT: tilt-in-space shower commode chair		40°
Kept in the wet room		Tilt-in-space
TILT-IN-SPACE	40°, with tilt-assist pedal	
SEAT HEIGHT	560 mm	
BACKREST	AdjustaBack, 50 cm, with flip-up armrests	
SEAT MODULE	Moulded, front opening	
FOOTPLATES	Adjustable 410–520 mm, folding plate	
CASTORS	Four dual-locking, 5 inch	
HEADREST	Moulded	
COMMUNE PAN	Autoclavable	
HIP BELT	Adjustable	

Mira Select Flex EV: shower valve		TMV3
		Thermostatic valve
TYPE	Thermostatic, TMV3	
FITTINGS	Long riser rail, 2 m hose and handset	
ACCREDITATION	RNIB Tried & Tested	

Drive DeVilbiss: shower stool		136 kg
		Safe working load

Geberit AquaClean Mera Care: wash-dry WC	
WASH	WhirlSpray shower, plus lady shower
DRY	Warm air dryer
CONTROLS	Remote control and touchless wall panel
ODOUR EXTRACTION	Fitted, with a soft-closing seat

Ropox Swing: height-adjustable washbasin		750–950 mm
		Basin height range
HEIGHT RANGE	750–950 mm, manual adjustment	
MAXIMUM USER WEIGHT	245 kg with the basin locked in the dock-in unit; 80 kg when it is swung out	
MOUNTING	Swing basin with dock-in unit	
TAP	Thermostatic mixer, lever operation, safe-touch body, TMV3 approved	

Kitchen

Worktops you can reach from a seated position. The living room has a rise-and-recline chair; rear French doors onto the patio give a 1720 mm opening.

NEFF Slide & Hide: oven	
Open-plan, wheelchair access to worktops	
OVEN	NEFF Slide & Hide: the door folds away underneath, so you are not reaching across a hot open door from a seated position
HOB	Gas, not induction: no electromagnetic field from the cooktop, which many guests with pacemakers prefer
STORAGE	Microwave and accessible storage at lower levels

What we can't promise

We can't guarantee every piece of specialist equipment at short notice. Some has to be hired in, and some depends on what's available that week. What we can promise is a straight answer quickly rather than leaving you hoping.

If you need a measurement that isn't here, ask and we will go and measure it.

Restwell Retreats, Whitstable, Kent. Optional home care comes from our sister company, Continuity of Care Services (CQC rated Good).

Phone [01622 809881](tel:01622809881) · Email hello@restwellretreats.co.uk · restwellretreats.co.uk/accessibility

This statement is a snapshot; the live accessibility page is always the current version.