

Access statement: equipment and measurements

A single-storey adapted bungalow in Whitstable, Kent, for disabled adults, families and carers. This is the spec sheet for occupational therapists, case managers and funding panels: makes, models, safe working loads and clearances, taken from the manufacturers' own spec sheets and our own tape measure.

965 mm Front door, clear opening	926 mm Internal doors, clear width	180 kg Ceiling track hoist, safe working load
--	--	---

IN THIS STATEMENT [Getting in](#) [Beds and hoists](#) [Wet room](#) [Kitchen](#) [What we can't promise](#)

Getting in

The driveway is resin-bound, level, and takes two cars including an adapted vehicle. There is no step up into the house from where you park.

Parking and getting around Single storey, no internal steps	
DRIVEWAY	Private, resin-bound: two off-road spaces, adapted vehicles welcome
STREET PARKING	Outside the house if you need extra room. No residents permit on this road; check signs on arrival in case street rules change
CAR TO FRONT DOOR	Step-free path
FLOORS	Level flooring throughout, no carpet lips
PATIO	Level hard-standing suitable for wheelchairs, with the enclosed garden and BBQ area beyond

Doorways and thresholds Clear openings, measured door-stop to door-stop		965 mm Front door clear opening
FRONT DOOR	965 mm clear, level threshold, no step	
PORCH OUTER DOORS	1720 mm opening, with a full-width ribbed entrance mat	
ALL INTERNAL DOORS	926 mm clear	
CONSERVATORY FRENCH DOORS	1720 mm opening	
REAR FRENCH DOORS	1720 mm opening, onto the patio	
PATIO THRESHOLD RAMP	20 mm high × 80 mm wide, non-slip	
DOOR FURNITURE	White doors with black lever handles, for contrast	
PORTABLE RAMPS	Fold-up ramps kept for the front door, yours for the stay and for days out	

Beds and hoists

The accessible bedroom can be set up with one or two profiling beds, depending on what you need.

Please bring your own slings

Slings are prescribed items of personal care. They are fitted to the individual and to the hoist mechanism, so we cannot safely supply them. Both hoists here take loop-style slings on a two-point spreader bar. Bring the sling the person already uses.

Amico GoLift 400: ceiling track hoist		180 kg
Main bedroom, over the profiling bed		Safe working load
SAFE WORKING LOAD	180 kg / 28 st	
LIFT AND LOWER	Powered	
TRAVERSE	Manual	
TROLLEY	Fixed	
ATTACHMENT	Lifting strap with loop	
COVERAGE	Full-room	
THOROUGH EXAMINATION	LOLER, every six months	

Accora CommunityBed		180 kg
Accessible bedroom, under the ceiling track		Maximum user weight
SAFE WORKING LOAD	215 kg / 34 st	
MAXIMUM USER WEIGHT	180 kg / 28 st	
MATTRESS PLATFORM	900 × 2000 mm	
OVERALL FRAME	930 × 2250 mm	
PLATFORM HEIGHT	220–800 mm	
PROFILING	Four-section: backrest with auto-regression, kneebreak, Trendelenburg and anti-Trendelenburg	
MATTRESS	Accora Allevia Comfort FirmEdge 90 cm, or Roma Medical MATT 1, both pressure-relieving	
SIDE RAILS	Fabric siderail set, 200 cm	
MINIMUM USER HEIGHT	1460 mm	
CARER ACCESS	Space for a carer on both sides of the bed	

Joerns Oxford Midi 180: mobile hoist		180 kg
Kept in the bungalow, usable in any room		Safe working load
SAFE WORKING LOAD	180 kg / 28 st	
TURNING RADIUS	1235 mm	
LEGS OPEN, EXTERNAL	1170 mm	
LEGS CLOSED, EXTERNAL	600 mm	
LEG HEIGHT	100 mm	

Joerns Oxford Midi 180: mobile hoist

Kept in the bungalow, usable in any room

180 kg

Safe working load

GROUND CLEARANCE	25 mm
SPREADER BAR HEIGHT	525–1660 mm
EMERGENCY DESCENT	Manual
THOROUGH EXAMINATION	LOLER, every six months. Last examined 14 May 2026

Bringing your own mobile hoist

The bed platform lowers to 220 mm. The Oxford Midi's legs stand 100 mm high with 25 mm of ground clearance, and its turning radius is 1235 mm. If you are bringing your own hoist, measure its leg height and turning circle against those figures, and tell us what you are bringing so we can set the room up for it.

AAL RS4: transfer and standing aid

185 kg

Safe working load

The RS4 is an Able Assist device. It is not an Arjo Sara Stedy, and the two are not interchangeable in a handling plan. Please check your care plan names the right one.

Wet room

Level-access throughout, with no lip at the entrance. Specified by Care Spaces.

Level-access wet room

Specified by Care Spaces

FLOOR	Level access with no lip, floor-level drain
TOILET SUPPORT RAIL	NYMAS hinged lift-and-lock, with drop-down leg
ASSISTANCE ALARM	Pull-cord with reset plate and over-door light
EXTRACTOR FAN	Fitted

RAZ-AT: tilt-in-space shower commode chair

Kept in the wet room

40°

Tilt-in-space

TILT-IN-SPACE	40°, with tilt-assist pedal
SEAT HEIGHT	560 mm
BACKREST	AdjustaBack, 50 cm, with flip-up armrests
SEAT MODULE	Moulded, front opening
FOOTPLATES	Adjustable 410–520 mm, folding plate
CASTORS	Four dual-locking, 5 inch
HEADREST	Moulded
COMMODE PAN	Autoclavable
HIP BELT	Adjustable

Mira Select Flex EV: shower valve

TMV3

Thermostatic valve

TYPE	Thermostatic, TMV3
FITTINGS	Long riser rail, 2 m hose and handset
ACCREDITATION	RNIB Tried & Tested

Drive DeVilbiss: shower stool

136 kg

Safe working load

Geberit AquaClean Mera Care: wash-dry WC

WASH	WhirlSpray shower, plus lady shower
DRY	Warm air dryer
CONTROLS	Remote control and touchless wall panel
ODOUR EXTRACTION	Fitted, with a soft-closing seat

Ropox Swing: height-adjustable washbasin

750–950 mm

Basin height range

HEIGHT RANGE	750–950 mm, manual adjustment
MOUNTING	Swing basin with dock-in unit
TAP	Thermostatic mixer, lever operation, safe-touch body, TMV3 approved

Kitchen

Worktops you can reach from a seated position. The living room has a rise-and-recline chair; rear French doors onto the patio give a 1720 mm opening.

NEFF Slide & Hide: oven

Open-plan, wheelchair access to worktops

OVEN	NEFF Slide & Hide: the door folds away underneath, so you are not reaching across a hot open door from a seated position
HOB	Gas, not induction: no electromagnetic field from the cooktop, which many guests with pacemakers prefer
STORAGE	Microwave and accessible storage at lower levels

What we can't promise

We can't guarantee every piece of specialist equipment at short notice. Some has to be hired in, and some depends on what's available that week. What we can promise is a straight answer quickly rather than leaving you hoping.

If you need a measurement that isn't here, ask and we will go and measure it.

Restwell Retreats, Whitstable, Kent. Optional home care comes from our sister company, Continuity of Care Services (CQC rated Good).

Phone [01622 809881](tel:01622 809881) · Email hello@restwellretreats.co.uk · restwellretreats.co.uk/accessibility

This statement is a snapshot; the live accessibility page is always the current version.